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# **BİYOMEDİKAL CİHAZ TEKNOLOJİLERİ**

**İNGİLİZCE SERVİS DOKÜMANLARI  
523EO0200**

**Ankara, 2011**

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INTRODUCTION	
Overview.....	v
Reference Documents.....	vi
Safety Information.....	vii
1.0 PRE-INSTALLATION.....	1-1
2.0 LUMISCAN INSTALLATION.....	2-1
3.0 SYSTEM POWER-UP AND OPERATION.....	3-1
4.0 THEORY OF OPERATION.....	4-1
5.0 SERVICE ADJUSTMENTS.....	5-1
6.0 PERIODIC MAINTENANCE.....	6-1
7.0 DIAGNOSTICS.....	7-1
8.0 REPLACEMENT PARTS.....	8-1
9.0 APPENDIX A.....	9-1
10. SCHEMATICS.....	10-1

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### 1.1.1.2. Kurulum Öncesi (Pre-Installation)

<b>5. PRE-INSTALLATION</b>	<b>24</b>
5.1 Mounting methods.....	25
5.2 Electric pre-setting .....	26
<b>6. INSTALLATION</b>	<b>27</b>
6.1 Timer set up (standard configuration).....	28
6.2 Assembling the mobile stand and timer installation .....	30
6.3 X-ray button .....	32
6.4 Assembling the arms .....	33
6.4.1 Assembling the extension arm.....	33
6.4.2 Assembling the scissors arm (DP arm).....	34
6.4.3 Assembling the stand arms set .....	36
6.5 Installation of the optional parts.....	37
6.5.1 External signaling devices.....	37
6.5.2 Installation of chemical screws.....	37
6.5.3 Counterplate .....	38
6.6 Electrical connection.....	39
6.6.1 Electrical connection for standard versions .....	40
6.6.2 Wall Support - Remote keyboard connection .....	42
6.7 Power-up sequence .....	43

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<b>5.2 Electric pre-setting</b>	
<ul style="list-style-type: none"><li>• Single-phase supply + ground</li><li>• Frequency</li><li>• Absorbed current during emission</li><li>• Apparent line resistance</li></ul>	<ul style="list-style-type: none"><li>230 V~ / 120 V~</li><li>50 Hz – 60 Hz</li><li>2.25 A / 4.5 A</li><li>0.8 Ω max / 0.4 Ω max.</li></ul>
	<b>NOTE:</b> <b>The device is intended for permanent installation.</b> <b>It is forbidden to connect the system by a plug in a wall socket for safety reasons.</b> <b>Mobile version has to be requested on order and it is not possible to transform wall mounted version into a mobile version in field.</b>
	<b>NOTE:</b> A circuit breaker with overcurrent protection must be connected to the intraoral X-ray equipment with the following features: <ul style="list-style-type: none"><li>- Nominal current: 6A</li><li>- Differential sensitivity: 0.03A.</li></ul>

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**Remove** the small box containing the rotor (Figure 2-1).

Remove the Power Cords and fuses.

Remove the top cardboard supports.


Remove the cardboard protector from the Cell Washer.

Grip the cell washer under the top edge of the front panel and under the top edge of the pump enclosure and lift it free of the carton. (Have someone hold the carton, if necessary).

**Place** the cell washer on a sturdy, level bench near a waste container or drain to receive decanted saline.

Carefully inspect the cell washer for external damage. Any damage discovered which is attributable to mishandling or shipping should be documented, and a signed inspection report should be furnished to the shipping company.

Complete the warranty registration card and return it to IEC.

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2	INSTALLATION .....	5
	2.1 Receiving Inspection .....	5
	2.2 Unpacking .....	5
	2.3 Power Configuration .....	7
	Voltage .....	7
	Fuses .....	7
	Power Cord .....	8
	2.4 Tubing Installation .....	8

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<b>2.2</b>	<b>Hardware Installation</b>
<b>2.2.1</b>	<b>Tools Required</b>
<b>2.2.2</b>	<b>AC Voltage</b>
<b>2.2.3</b>	<b>Power Cable</b>
<b>2.2.4</b>	<b>Media Support Assembly</b>
<b>2.2.5</b>	<b>Screen Exit Assembly</b>
<b>2.2.6</b>	<b>Installation of Data Control Board (DCB)</b>
<b>2.2.7</b>	<b>Installation of Control Interconnect Cable</b>
<b>2.2.8</b>	<b>Installation of the ACR-2000 Eraser</b>

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### 2.2.3 Power Cable

The LUMISCAN ACR-2000 READER utilizes an international IEC grade connector for the power cable. Systems are shipped with a standard NEMA 5-15 hospital grade cable. This cable requires replacement depending upon the country of installation. Insert the female end into the input power socket.

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### 2.3 Driver Installation

#### 2.3.1 Corection LUT Files

#### 2.3.2 DOS, Windows 3x and Windows 95 and windows 98 Driver Installation

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<b>Chapter 4: Setup and Operation</b>	
4.1 System Controls .....	21
4.2 System Components .....	22
4.3 Setup .....	23
4.4 Touchpad .....	24
4.5 Accessories .....	24
<b>4.6 Preparing the System for Operation .....</b>	<b>25</b>
4.6.1 Installing and Removing the Battery .....	25
4.6.2 Using AC Power/Charging Battery .....	26
4.6.3 Connecting to AC Power .....	27
4.6.4 Connecting and Removing Transducers .....	28
4.6.5 Turning the System On and Off .....	28

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<b>2.0 Installation</b>	<b>12</b>
2.1 Probe mounting	12
2.2 Opening the analyzer’s enclosure (wall-mount version)	13
2.3 Mounting the analyzer (wall-mount version)	13
2.4 Wiring	14
2.5 Connecting the probe	16
2.6 Power wiring	17
2.6.1 AC supply	17
2.6.2 DC supply	17
2.7 Connecting output signals	18
2.7.1 Relays	18
2.7.2 Current outputs	18
2.7.3 Voltage outputs	18
2.8 Auxiliary input signals	19
2.8.1 4 to 20 mA loop powered device	19
2.8.2 Active 4 to 20 mA device	19
<b>3.0 Programming and Operation</b>	<b>20</b>
<b>3.1 Applying power</b>	<b>20</b>

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### 3.1 Applying power

When all wiring has been secured, power may be applied to the analyzer. (No power switch is supplied on the 2350—the system turns on as soon as power is applied.) Upon **power-up** the analyser performs a software diagnostics test, and then displays humidity and temperature values.

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PERIODIC MAINTENANCE .....	34
<b>Preventative Maintenance</b> .....	34
Cleaning Instruction: .....	34
Cleaning instruction .....	35

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### ***Preventative Maintenance***

Preventive maintenance inspections should be performed on a semi-annual basis and should include the following:

- Inspect all electrical wires and cables for abrasions, cuts or cracks.
- Inspect all control buttons and knobs to ensure that they are secured properly and accurately.
- Inspect all labels and indicators on the operator's control panel to ensure legibility and accurate alignment.
- Confirm mA, kVp and timer accuracy and adjust as necessary.

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Section 8: Routine Maintenance and Shipment.....	41
8.1 Routine Check of Backup Systems .....	41
8.2 Hand Pump Test.....	41
8.3 Self-Test Routine .....	41
8.4 Emergency System Test .....	41
8.5 Packing for Shipment .....	41
8.6 – <u>For BVS® Blood Pumps and AB5000™ Ventricles</u> – Cleaning .....	41
8.7 Scheduled Maintenance .....	42
8.8 Replacement Parts .....	42
8.9 Service.....	42
8.10 _Ordering Information.....	43

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Maintenance	Period
Perform cleaning, functional test, and calibration.	1 year
Replace batteries. *	2 years
Replace solenoid valves.	5,000 hours
Replace compressor.	5,000 hours
Replace proportional valves.	5,000 hours
Estimated useful Console life. **	10,000 hours

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SERVICE ADJUSTMENTS .....	29
SCHEMATICS .....	30
DIAGNOSTICS .....	31
Troubleshooting SCR Block Gates .....	31
Troubleshooting Main SCR Gate Drive .....	31
Troubleshooting the Backup SCR Gate Drive .....	31
Renewal Parts .....	33

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Service Provider					
Name:					Date:
Company:					
Address:					
Phone Number:			Fax Number:		
E-mail address:					
Device Description					
Name:			Serial Number:		
Part Number:		Lot Number:		Revision:	
Software Version:			Other Identifiers:		
Event Description					
Diagnosis					
Service Performed					
Performed By:					Date:
Actions:					
Parts Removed					
Part Name	Part Number	Serial Number	Lot Number	Rev	Replaced By
Parts Installed					
Part Name	Part Number	Serial Number	Lot Number	Rev	Replaced By
Tests Performed (attach test data)					
Test:			Test:		
Performed By:			Performed By:		
Result: Pass <input type="checkbox"/> Fail <input type="checkbox"/>			Result: Pass <input type="checkbox"/> Fail <input type="checkbox"/>		

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<b>IV. Care and Maintenance</b>	
STIMUPLEX Dig RC Short Test .....	.16
To Check Remote Control .....	.17
Cleaning .....	.17
Battery replacement .....	.17
Replacement of the Electrode Cable .....	.18
Instructions for Electrode Cable Replacement .....	.18
<b>Troubleshooting</b> .....	.19
Service and Technical Support .....	.20

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<b>Troubleshooting:</b>	
<b>PROBLEM:</b>	Yellow light (LED) does not light during short test.
<b>SOLUTION:</b>	See technical testing on page 21 of this manual.
<b>PROBLEM:</b>	During the procedure, the digital display flashes and the signal tone does not change its length.
<b>SOLUTION:</b>	(i) Check skin electrode to confirm you have sufficient contact with the patient's skin. (ii) Check that the alligator clip is secured to the skin electrode (iii) Check that there is a secure connection between the needle and the cable. Perform the technical testing (see page 21) to confirm cable operation.
<b>PROBLEM:</b>	Red "Bat" LED flashes (continuously).
<b>SOLUTION:</b>	The red LED will flash if battery voltage drops below 6 volts. A blockade which is in progress can be continued, but the battery must be replaced afterwards.
<b>PROBLEM:</b>	There is no audio signal after the unit is switched on.
<b>SOLUTION:</b>	The battery is low, replace the battery.

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8.	TRUBLE-SHOOTING GUIDE .....	8.1
	Basic Problems .....	8.1
	Monitor .....	8.1
	Stimulator .....	8.3
	Patient 8.4	
	WR Address and Phone/FAX Numbers .....	8.4

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<b>1. SYMPTOM: The Stimulator does not work</b>
<input type="checkbox"/> <i>Is the power on (indicated by a lighted switch)?</i>
NO: Turn the power on by pressing the POWER SWITCH. The switch should light.
<input type="checkbox"/> <i>Is the REPLACE BATTERY light on?</i>
YES: Replace the battery pack with a freshly charged battery pack. Make sure that the battery pack that you have removed is put in the battery charger to charge.
<input type="checkbox"/> <i>Is the PULSE light blinking? Is the CURRENT INTENSITY display lighted?</i>
NO: Stimulator may need repair. Contact the factory for advice.

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Hex Code	Explanation	Recommended Action
A	Measured value of isolated reference supply on front end was low	<ol style="list-style-type: none"> <li>1. Check power supply.</li> <li>2. Replace main PCB.</li> </ol>
B	On a NPB-4000 rev 4 PCB or earlier. The measured value of isolated reference supply on front end was high.	<ol style="list-style-type: none"> <li>1. Check power supply.</li> <li>2. Replace main PCB.</li> </ol>
C	On a NPB-4000 rev 5 or later, or on NPB-4000C. The measured value of isolated reference supply on front end was high.	<ol style="list-style-type: none"> <li>1. Check power supply.</li> <li>2. Replace main PCB.</li> </ol>
64	MP-205 SpO2 module is not sending messages to host CPU.	<ol style="list-style-type: none"> <li>1. Ensure MP-205 module is properly connected.</li> <li>2. Replace MP-205 module.</li> <li>3. If problem persists, replace main PCB.</li> </ol>
65, 66	MP-205 SpO2 module has detected an error during initialization.	<ol style="list-style-type: none"> <li>1. Replace MP-205 module.</li> <li>2. If problem persists, replace main PCB.</li> </ol>
6E-71	MP-205 SpO2 module has detected an error on its serial port.	<ol style="list-style-type: none"> <li>1. Replace MP-205 module.</li> <li>2. If problem persists, replace main PCB.</li> </ol>

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
<b>Display Replacement</b>	
1	Set the new display in place.
2	Install the four screws that hold the Display in place. Torque the screws to 5.5 inch pounds.
3	Connect the two connectors that connect the Display to the Main PCBA.
4	Place the Control Panel in place.
5	Reinstall the four screws that hold the Control Panel in place. Torque the screws to 5.5 inch pounds.

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
SERVICE ADJUSTMENTS.....	29
<b>SCHEMATICS.....</b>	<b>30</b>
DIAGNOSTICS .....	31
Troubleshooting SCR Block Gates.....	31
Troubleshooting Main SCR Gate Drive .....	31
Troubleshooting the Backup SCR Gate Drive.....	31
Renewal Parts .....	33

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<h3>5.0 Calibration of the Interface Unit</h3> <p>The calibration of the EPIC-4000 interface consists of two procedures:</p> <ol style="list-style-type: none"> <li>1. Calibration of gain for each channel</li> <li>2. Calibration of offsets for each channel</li> </ol> <p><b>5.1 Patient Amplifier</b></p> <p>Each Patient Amplifier requires adjustment of the Driven Ground DC Offset. One trim pot only (R51) adjusts the offsets for both channels. Follow this procedure for each Patient Amplifier board in the system:</p> <ol style="list-style-type: none"> <li>1) Place shorting jumpers across the positive (+) and negative (-) inputs for channels 1 and 2 on the patient cable.</li> <li>2) Connect the negative (-) lead of an oscilloscope on pin 2 of the Burr-Brown 722 DC/DC Converter IC (Isolated Ground).</li> <li>3) Connect the oscilloscope test probe on the Patient Cable Ground (<b>Gnd</b>) lead,</li> <li>4) Make sure that no 50 or 60 Hz interferences are present. Then adjust the trim pot R51 so that at the most sensitive setting of the oscilloscope (2 mV / div is sufficient), the DC offset is adjusted to 0 mV.</li> </ol>
--

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#### 2.4 SG-2002 Ganzfeld Stimulator

The SG-2002 Ganzfeld Stimulator consists of two units: the control box and Ganzfeld. The SG-2002 is controlled by the system computer via a serial port connection.

Xenon flash tubes are used to produce the required stimulus. Any single flash or flicker can be triggered with or without background light. The intensity of flash or flicker can be set to any of three dB levels, +25dB, 0dB or -25dB. The brightness of the background light is pre-set to any of three levels for flash/flicker or EOG. Three red LEDs, either set to bright or dim brightness, are mounted at the back of the Ganzfeld for EOG testing.

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#### II. Product Description

Design and Functions .....	7
Remote Control .....	9
Technical Data .....	10

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
OPERATION	7
Control Panel Controls & Functions ....	7
Operating Procedure .....	8

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1. **ON Button** - Pushing this button will turn the unit ON.
2. **OFF Button** - Pushing this button will turn the unit OFF.
3. **BUCKY Switch** - Selecting Table or Wall will allow the respective optional Bucky to energize when an exposure is initiated. Selecting Off will inhibit either Bucky from energizing when an exposure is initiated.
4. **Minor kVp** - This rotary switch selects kVp in 1-2 kVp increments and is used in conjunction with the Major kVp switch to select the desired kVp.
5. **Major kVp** - This rotary switch selects kVp in approximately 10 kVp increments and is used in conjunction with the Minor kVp switch to select the desired kVp.
6. **mA Selector** - This rotary switch selects the desired mA level with the corresponding small or large focal spot. S after the number indicates small focal spot and L indicates large focal spot.
7. **Time Selector** - This rotary switch selects the desired timer setting.
8. **PREP Button** - Pushing and holding this button will initiate the Prepare cycle. The anode in the x-ray tube begins to spin and the filament in the

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<b>Chapter 2: Safety</b>		
2.1	Electrical Safety .....	5
2.2	Equipment Safety .....	6
2.3	Battery Safety .....	6
2.4	Biological Safety .....	7
2.5	Labeling Symbols .....	7

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## SAFETY

Please read this information carefully prior to installing or using this equipment.

1. The unit described in this manual is designed to be operated only by trained personnel. Any adjustments, maintenance or repair must be carried out as defined in this manual, by a person qualified to be aware of the hazards involved.
2. It is essential that both operating and service personnel employ a safe system of work, in addition to the detailed instructions specified in this manual.
3. The covers of the unit should only be removed by personnel who have been trained to avoid the risk of shock.
4. References should always be made to the Health & Safety data supplied with any chemicals used. Generally accepted laboratory procedures for the safe handling of chemicals should be employed.
5. If it is suspected that safety protection has been impaired in any way, the unit must be made inoperative and secured against any intended operation. The fault condition should immediately be reported to the appropriate servicing authority.

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<b>SAFETY ASPECTS</b>	
2.1	Warnings .....
2.2	Protection from X-rays .....
2.3	<b>Environmental risks and disposal.....</b>
2.4	Symbols in use .....
2.5	Attention points during installation ..

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## Environmental risks and disposal

Some parts of the equipment contain material and fluids which must be disposed of in special areas designated by the local health authorities at the end of the equipment's life cycle.

In particular the equipment contains the following materials and / or components:

- **Tubehead:** external packages in non-biodegradable plastic, dielectric oil, lead, copper, brass, aluminium, resin, tungsten, beryllium
- **Power supply and remote control:** external packages in non biodegradable plastic, iron, copper, plastic reinforced by fibre glass
- **Tubehead extension:** iron, aluminium, copper.

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<b>10. SERVICE ACTIONS REQUIRING PART REPLACEMENTS</b>	<b>65</b>
10.1 Replacement of the tubehead.....	65
10.2 Replacing the generator card .....	67
10.3 Replacing the logic card .....	68
10.4 Replacing the keyboard.....	69
10.5 Replacing the scissors arm .....	70
10.5.1 Replacing the scissors arm with wall Timer (standard configuration)...	70
10.5.2 Replacing the scissors arm with Timer with mobile stand .....	71
<b>11. MAINTENANCE</b>	<b>72</b>
<b>12. SCHEMATICS AND DRAWINGS</b>	<b>73</b>
<b>13. SPARE PARTS</b>	<b>81</b>

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
<b>1</b>	<b>General Notes</b>
	General notes . . . . .
	Safety . . . . .
	System configuration . . . . .
<b>2</b>	<b>Room Planning</b>
	Dimensions: MOBILETT Plus HP . . . . .
<b>3</b>	<b>Technical Data</b>
	Electrical data . . . . .
	Weight and heat dissipation . . . . .
	<b>Environmental conditions</b> . . . . .
	Packing and transport routes . . . . .
	Surface colors . . . . .
<b>4</b>	<b>Changes to Previous Version</b>

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Chapter 4: Setup and Operation	
4.1 System Controls .....	21
4.2 System Components .....	22
4.3 Setup .....	23
4.4 Touchpad .....	24
4.5 Accessories .....	24
4.6 Preparing the System for Operation .....	25
4.6.1 Installing and Removing the Battery .....	25
4.6.2 Using AC Power/Charging Battery .....	26
4.6.3 Connecting to AC Power .....	27
4.6.4 Connecting and Removing Transducers .....	28
4.6.5 Turning the System On and Off .....	28
<b>4.7 Upgrading the System Software .....</b>	<b>29</b>
4.7.1 Obtaining a License Key .....	34
4.7.2 Installing a License Key .....	34
4.7.3 To Display the System Information Screen .....	35
4.7.4 To Display the License Update Screen .....	35

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**To upgrade the system software:**

**Caution:** To prevent loss of data or loss of images, transfer all images on the CompactFlash card before performing the upgrade. Images remaining on the CompactFlash card cannot be viewed or deleted after the upgrade.

*Note: If you use SiteLink Image Manager for transferring images, SiteLink 2.1 is required after the upgrade. Contact the Technical Support Department to receive a new version of the program.*

- 1 Remove any transducer or Triple Transducer Connect from the Titan system.
- 2 Connect the Titan system directly to the power supply or through the docking system/mini-dock. See “Connectivity” on page 73.
- 3 Insert the CompactFlash card into the back slot.  
The system displays the following message:

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## 9. Glossary

**contralateral shoulder:** The shoulder of the patient on the opposite side of the site of surgery. This is the location of the *stimulator return electrode* when a *monopolar stimulator probe* is used.

**current intensity:** The amount of current being delivered by the *stimulator*.

**electrosurgical detection circuit:** A circuit in the Brackmann II that detects the radio signals characteristic of an electrosurgical device, such as electrocautery. When the signal is detected, the speaker of the Brackmann II is muted, eliminating the loud static (*artifacts*) that would otherwise occur during electrocautery. The sensitivity is set by the MUTE SENSITIVITY dial.

**EMG signal (activity):** Electromyographic signal (activity). Electrical energy resulting from muscle movement, or minute electrical signals or activity in patient tissue that is amplified and then displayed and/or broadcast on an EMG device such as the Brackmann II.

**ELECTRODE WARNING indicator:** An indicator that lights when there is improper *impedance* in the *bipolar, hooked-wire electrodes* used to monitor *EMG activity*.

**electrode:** A fixed electrical connection used for conducting electrical current from patient tissue.

**impedance:** As relating to the Brackmann II, obstruction in the flow of electric current through patient tissue.

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<b>Index</b>	
<b>T</b>	<b>V</b>
technical support, contact 67	VCR problem 39
THI 21	video output tests 66
touchpad 22	video recording 22
transducer	
acquisition subsystem 10	
problems 40	
specifications 16, 18	
storage and shipping 18	
troubleshooting symptoms table 39	
	<b>W</b>
<b>U</b>	warnings 5
update 22	
	<b>Z</b>
	zoom
	description 21
	problem 39

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### 5.1 Replacing the mains fuse on the microscope

#### **CAUTION**

ALWAYS DISCONNECT YOUR EQUIPMENT BEFORE DOING ANY REPAIR.

#### Location of Mains Fuse

The mains fuse of the TC Series Inverted Microscopes is located on the back panel of the microscope as shown.

#### Instructions to replace the mains fuse:

- Turn the power switch to the off position.
- Unplug the microscope from the wall outlet.
- Loosen the fuse cap from the fuse holder.
- Remove the blown fuse from the fuse cap.
- Replace the fuse with the CORRECT type and rating which is: IEC Standard 5 X 20mm 3amp Fast-Acting such as Littelfuse 217 Series or Buszmann GDB Series
- Reinstall the fuse holder with the new fuse installed.

#### **DANGER**

NEVER USE REPLACEMENT FUSES OF A DIFFERENT RATING.

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<b>2. SAFETY ASPECTS</b>	<b>2</b>
2.1 Warnings .....	3
2.2 Protection from X-rays .....	4
2.3 Environmental risks and disposal.....	5
2.4 Symbols in use .....	6
2.5 Attention points during installation .....	7

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6.4	Assembling the arms .....	33
6.4.1	Assembling the extension arm.....	33
6.4.2	Assembling the scissors arm (DP arm).....	34
6.4.3	Assembling the stand arms set .....	36
6.5	Installation of the optional parts .....	37
6.5.1	External signaling devices.....	37
6.5.2	Installation of chemical screws.....	37
6.5.3	Counterplate .....	38
6.6	Electrical connection.....	39
6.6.1	Electrical connection for standard versions .....	40
6.6.2	Wall Support - Remote keyboard connection .....	42
6.7	Power-up sequence .....	43

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At start-up the timer automatically checks all the internal hardware and in particular:

- Display: all segments light for 2 seconds
- LED: they all light for 2 seconds
- CHECKSUM check of the storage program and display of the words "CHS" if the check result is negative or in the event of non-initialized memory
- Buzzer check (rings)
- Check of all buttons and X-ray buttons in particular, no button must be pressed. If any buttons are pressed one of the LEDs corresponding to this key will light up and an error message will appear
- Display of software version (ex. 1.01).

- A) Electrical connection
- B) Power up sequence
- C) Assembling the arms
- D) Installation of the optional parts
- E) Hiçbiri

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Check the red and green stimulator battery indicators for the battery condition. A green indication indicates adequate power remains. A red indicator indicates that the stimulator batteries must be replaced.

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The Remote Probe, Remote Indicator Light, and the Muscle Movement Sensor must be gas sterilized only.

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*The use of adapters to plug in more than one instrument to the same outlet (power strip) does not provide a reliable power supply. It may produce overheating and create a potential fire hazard. ECG electrical instrument directives forbid such practices.*

3. Turn the power on.
4. Select the preferred Language on the initial screen.
5. From the **Main** screen select **Service/ Edit Setup**. Set the elevation to the appropriate value for the location of the laboratory. (See Chapter 12 for more details.)

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### **Equipment Safety**

To protect your ultrasound system, transducer, and accessories, follow these precautions.

**Caution:** Excessive bending or twisting of cables can cause a failure or intermittent operation.

To avoid damaging the power supply, verify the power supply input is within the correct voltage range. See “[Electrical](#)” on page 18 in [Chapter 3](#).

Improper cleaning or disinfecting of any part of the system can cause permanent damage. For cleaning and disinfecting instructions, see [Chapter 5, “Cleaning and Disinfecting.”](#)

Do not use solvents such as thinner or benzene, or abrasive cleaners on any part of the system.

Remove the battery from the system if the system is not likely to be used for some time.

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


Procedure	
<i>To remove power supply assembly:</i>	
a) On rear panel of monitor, remove eight screws fastening power supply assembly.	
Illustration	
 <p>Screws fastening power supply assembly</p>	
b) Carefully lift power supply assembly from rear case.	
c) Power supply assembly may be disassembled into the following elements:	
<ul style="list-style-type: none"> <li>• Power supply PCB</li> <li>• NIBP pump</li> <li>• Heat sink/chassis</li> <li>• Main ribbon cable</li> </ul>	
Illustration	
	

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2.0 MAINTENANCE SCHEDULE	
Detailed maintenance procedures are in sections 3 and 4.	
<b>2.1 THE FOLLOWING ACTIONS MUST BE TAKEN AFTER EACH USE:</b>	
1. Sterilize equipment.	emir cümlesi (imperative)
2. Recharge monitor batteries	
3. Check stimulator batteries and replace if necessary.	
4. Check front panel controls.	
<b>2.2 THE FOLLOWING ACTIONS MUST BE TAKEN EVERY THREE MONTHS:</b>	
1. Check monitor calibration.	
2. Check stimulator calibration.	
3. Replace stimulator batteries.	
4. Clean the exterior with a water damp cloth.	
<b>2.3 THE FOLLOWING ACTIONS MUST BE TAKEN EVERY YEAR:</b>	
1. Unit should be returned to WR Medical Electronics every year for routine calibration.	
2. Replace Nickel cadmium monitor batteries.	

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<b>Symptom</b>	<b>Solution</b>
<u>System will not power on.</u>	Check all power connections. ✓ Perform the following sequence: remove DC input connector and battery; wait 10 seconds; connect DC input or install battery; press the power key. ✓ Ensure the battery is charged. ✓
System image quality is poor.	Adjust the LCD screen to improve viewing angle. Adjust the brightness, as necessary, to improve image quality. Adjust the gain.
Zoom does not work.	Press Freeze key. Zoom does not work when the image is frozen.
No CPD image.	Adjust the gain.
No DCPD image.	Adjust the gain.

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### **Assembling the mobile stand and timer installation**

1. Cross the two base tubes (1) into the provided cut, fixing them together with the screw (2) and relevant nut (3). Do not tighten the screw completely.
2. Position the base plate (4) and fix it with the four screws (5). If necessary, reposition the two base tubes (1) slightly in order to align the relevant holes on the plate.
3. Lock the nut (3) in order to block the base tubes (1) permanently.
4. Assemble the stand column (6) on the base plate (4) with the four screws (7).


Figure 6-2

5. Assemble the timer support plate (8) to the stand column, fixing it with the two screws (9).


#### **WARNING:**

The timer must be fixed to the support plate after assembling the scissors arm (see paragraph 6.4.3).

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TESTING FOR STIMULUS OUTPUT: STIMULUS VERIFICATION
The unit <b>may be</b> tested for output by touching the active electrode to the ground pad. If the amber INCOMPLETE STIMULATION light goes out, the instrument is working properly.
The stimulus verification circuitry monitors the integrity of the patient/instrument circuit, including continuity of the probe and cable. When proper current is being delivered, the amber INCOMPLETE STIMULATION light will go out. It will flash at all other times - i.e. broken cables, disconnected cables, poor ground, poor probe contact, etc.
Verification of stimulation <b>can also be</b> obtained by touching the forearm or wrist area of a test subject, starting with a low current setting and increasing to a reasonable level (up to 6 or 8 mA may be required due to the narrow pulse width of .0002 sec.). Use a non-sterile probe for this procedure so that sterile probes will be available for surgery.
The unit <b>may also be</b> tested on an oscilloscope using a 1K precision resistor across the output. The oscilloscope will then display the pulse amplitude directly in milliamperes (1 volt = 1 mA).
The active probe <b>can be</b> shorted to the indifferent electrode without damage to the circuitry.

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CHECKING THE SENSOR CIRCUIT
You <b>may test</b> the functioning of the monitoring circuit and audible signal by lightly touching the sensor. When the sensitivity switch is set to “1”, the sensor will be so sensitive that if you set it on a flat surface and blow on it, you will set it off. You should be able to touch the patient’s cheek lightly and get a response.
CARE OF THE SENSOR
The sensor is quite delicate and <b>should be treated</b> with care. Pulling on wires, repeated bending, especially sharp bends, can cause broken wires or intermittent false signals. The sensor does not need to be sterilized because the patient’s mouth is not a sterile field. The sensor may be cleaned by wiping with cold cleaning or sterilizing solutions and may be gas autoclaved. It must not be immersed in cleaning solutions and must not be steam autoclaved because such procedures will draw moisture into the sensor and the cable, causing corrosion and malfunctioning of the instrument.

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The equipment gives optimum performance with minimal maintenance.  
Subject Verb(s) Object

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1. “Proportional Valve Display” indicates whether the proportional valve is open or closed.

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2. The equipment requires no routine service or calibration other than cleaning and battery maintenance.

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3. The lamp and lamp house become extremely hot during and after operation.

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4. This manual doesn't contain information concerning the installation, operation and maintenance of biological microscopes.

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5. The unit requires a grounded power supply.

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2. The equipment is designed to give optimum performance with minimal maintenance.

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3. The Power Switch is located at the front lower-right of the cabinet.

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4. The unit is provided with two power cords.

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1. If such an error occurs during monitoring operation, an audible alarm tone sounds.

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2. If saline flow is interrupted or depleted, the run stops, an alarm sounds and the ADD SALINE LED lights.

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1. Clean all external surfaces with this solution.

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2. Press the **ENTER** button to clear all alarms.

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3. Perform an internal test by pressing the **TEST** button.

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4. Remove all the screws located on the rear of the console base.

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### 2.5.3.3. Koşullu Emir Cümleleri (Imperatives with Conditionals)

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If + (genellikle present ya da bazen past tense), imperative

Replace the inspect filter if damaged.

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If the information in this manual is insufficient, contact us.

if clause (if + present tense) imperative

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1. If you can not read any impedances, verify that electrodes are connected properly.

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2. If you do not see a calibration pulse after reconnecting the amplifier cables, verify that the computer is **set** up correctly.

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3. If all stimulators are functioning properly, continue with the hard disk verification.

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4. If all amplifier channels carry noise, continue with step c.

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Verb (V) + complement to do something

Buna göre;

To stop pumping, press the appropriate OFF button twice within 13 seconds.

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1. To replace the motor unit, remove the cable between 2-pin terminals and motor connector.

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To reduce the vacuum level, press the function (Func) key.

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2. To return to normal operation, press and hold the ACTIVATE button.

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Unscrew the screw clamp enough to fit the clamp around the metal post.

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1. Always disconnect the power cord from the back of the microscope when not being used, when cleaning your instrument or when making any repairs. (When bağlacıyla kurulan emir cümlesi)

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Press the up arrow until the target flow setting stops increasing.

Hold the Bed Mount so the metal post faces up.

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**It is recommended that.....(present tense) (.....tavsiye edilir.)**

**We recommend that.....(present tense) (.....tavsiye ederiz.)**

**You are recommended that....(present tense)**

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1. It is recommended that you use it with safety glasses.

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2. We recommend that the instrument be returned to the factory for servicing.

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We advise that....(present tense)

You are advised that....(present tense)

It is advisable to use a grounding cable G/V with section 2.5 mm<sup>2</sup>.

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3. You are advised not to switch the equipment off immediately after an exposure.

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1. The equipment should be cleaned daily.

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2. You must be aware of all safety issues when you install and operate the equipment.

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3. Testing of the alarm should be done at least every six months

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4. The on/off key must be depressed for 3 seconds to turn off the device.

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The equipment may be powered by an external battery.  
Subject                      verb3

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1. This can damage the equipment or other devices.

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2. The delay period can be set from 0 to 30 minutes.

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3. The device can be operated using either AC or DC power.

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4. The exterior surfaces can be cleaned with a soft cloth.

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5. A simple outlet tester can be used to correct outlet polarity and the presence of a grounded circuit.

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6. This could possibly lead to the instrument no longer performing accurately.

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During the Weaning operation, the Console's displayed flow may vary for several beats from the selected flow setting.

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#### **Örnek 1:**

If the CIRCUIT FAULT alarm sounds immediately after the ENTER button is pressed:

- The humidifier door may not be closed properly.
- The cartridge/circuit may be faulty.

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### Örnek 2:

**Never disconnect or reconnect  
the amplifier or amplifier arm  
with the system power on.  
This practice will damage the  
amplifiers.**

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
1. This has occurred when patient conditions change and/or when the Console periodically adjusts Blood Pump ejection duration to optimize flow.
2. As long as the Console line cord is plugged into an AC power source, the internal battery will be kept in a charged state.

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<ol style="list-style-type: none"> <li>1. Insert the Humidifier Cartridge into the cartridge door. Latch the door shut to make all electrical connections between the cartridge and the Life Pulse</li> <li>2. Connect the green gas inlet tube to the large barbed connector labeled GASOUT.</li> <li>3. Connect the small diameter purge tube to the barbed connector labeled PURGE.</li> <li>4. Locate the water inlet tube, an opaque tube running from the lower right corner of the cartridge with a check valve and a Leur fitting on the end.</li> <li>5. Install the water inlet tube by placing the opaque tube into the pump housing, closing the pump door securely to pinch the water inlet tube inside, and connect the water transfer tube to the Luer connector on the end of the water inlet tube.</li> <li>6. Tap the sterile water bag with the other end of the water transfer tube and open the tubing clamp.</li> </ol> <p>Use only sterile water for the cartridge water supply.</p>	
---	--

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
Weight of equipment and detachable parts	
Gross weight including packing	35 kg
Net weight of equipment in standard configuration	22 kg
60 cm extension arm (standard)	2.9 kg
80 cm extension arm	3.5 kg
30 cm extension arm	1.9 kg
Scissors arm	9 kg
Wall plate with generator	5 kg
Tubehead	4.5 kg

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Hex Code	Explanation	Recommended Action
A	Measured value of isolated reference supply on front end was low	1. Check power supply. 2. Replace main PCB.
B	On a NPB-4000 rev 4 PCB or earlier. The measured value of isolated reference supply on front end was high.	1. Check power supply. 2. Replace main PCB.
C	On a NPB-4000 rev 5 or later, or on NPB-4000C. The measured value of isolated reference supply on front end was high.	1. Check power supply. 2. Replace main PCB.
64	MP-205 SpO2 module is not sending messages to host CPU.	1. Ensure MP-205 module is properly connected. 2. Replace MP-205 module. 3. If problem persists, replace main PCB.
65, 66	MP-205 SpO2 module has detected an error during initialization.	1. Replace MP-205 module. 2. If problem persists, replace main PCB.
6E-71	MP-205 SpO2 module has detected an error on its serial port.	1. Replace MP-205 module. 2. If problem persists, replace main PCB.

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- 1 Remove the four black pads located on bottom of unit.
- 2 Unscrew the four screws of the bottom panel with a screwdriver for recessed-head screws.
- 3 Pull off the bottom panel carefully upwards and position it by the side of the device (be careful of the battery connecting cable).
- 4 Disconnect the electrode cable plug from the plate and take the cable out of the clamping arrangement (clamped joint) of the front panel.
- 5 Press the new electrode cable with its rubber funnel (rubber bushing) into the clamping arrangement of the front panel and adjust the cable position.
- 6 Attach the plug of the electrode cable to the counter plug/socket on the plate. Don't use excessive force and be careful of the codings.
- 7 Replace the bottom panel carefully to seal. Take care that no cables are being compressed (battery cable).
- 8 Screw in all four screws.
- 9 Carry out an operation check.

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**CAUTION: The following safety precautions should be observed prior to attempting the light source replacement procedure.**

- 1. DISCONNECT THE UNIT FROM THE MAINS SUPPLY PRIOR TO REMOVING THE TOP COVER.**
- 2. SAFETY GLASSES MUST BE WORN WHEN UV EMISSIONS ARE PRESENT.**
- 3. DO NOT TOUCH THE LIGHT OUTPUT MODULE OF THE LAMP MODULE.**

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Condition	Recommended Action
1. The NPB-4000/C fails to power-up when the On/Standby switch is pressed.	<ol style="list-style-type: none"><li>1. Ensure power cord is plugged into operational AC outlet of appropriate voltage and frequency. Ensure green BATTERY CHARGING/AC SOURCE indicator is lit. If indicator is not lit, replace power supply assembly.</li><li>2. Check fuses located on power supply assembly above AC inlet receptacle. Replace fuses if necessary.</li><li>3. Inside monitor, check main ribbon cable and ensure that it is connected to main PCB and power supply assembly.</li><li>4. Ensure keypad is plugged into Main PCB. If connection is good, replace keypad.</li><li>5. If the problem persists, replace main PCB.</li></ol>

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